



KN95 Disposable Facepiece

Particulate Respirator

20pcs/box

GB2626-2006 KN95
Nelson Lab Reports Available



Product Name: KN95 Disposable Facepiece

Main Material: Spunbond Non-woven and Meltblown Filter Fiber

Product Standard: GB2626-2006 KN95

Specifications: DM9501

Packaging: 20pcs/box

Life Time: 3 Years (Please see the indications on the package for production batch number and expiration date)



Email: export@imdmedical.com
Find us on Facebook and Instagram with
IMDMEDICAL




Innovative Material and Devices, Inc.

Address: 3F, Building #5, No. 615 Fengdeng Rd, Jiading District, Shanghai 201801, China

Toll free: +86 400 8202 169 Web: www.imdmedical.com

– 4 Ply Guarantee · Upgraded Protection –

Filtering Efficiency \geq 95%



High Quality Non-woven Spunbond Fiber

Effective isolation away from harmful substances

Thickened Electrostatic Absorption Melt-blown Filtering Fiber

A powerful protection against both oily and non-oily particles

Ultrafine Fiber Melt-blown Filtering Layer

Low respiratory resistance

Soft Non-woven Spunbond Fiber

Moisture resistant and comfortable to wear

Breathing Healthy



Droplet Protection



Anti-Dust



Comfort



Filter Efficiency \geq



Innovative Material and Devices, Inc.


Address: 3F, Building #5, No. 615 Fengdeng Rd, Jiading District, Shanghai 201801, China

Toll free: +86 400 8202 169 Web: www.imdmedical.com

– Effective Filtration · Smooth Breathing –

Thickened Electrostatic Absorption Melt-blown Filter Fiber

Two layers high quality melt-blown filter fiber maintain at least 95% filtering efficiency and low respiratory resistance



Exhalation Resistance(160L/min)

≤ **250** Pa

Inhalation Resistance(95L/min)

≤ **350** Pa

Filtering Efficiency

≥ **95** %

Total Inward Leakage

11 %



Innovative Material and Devices, Inc.

Address: 3F, Building #5, No. 615 Fengdeng Rd, Jiading District, Shanghai 201801, China

Toll free: +86 400 8202 169 Web: www.imdmedical.com

– Ergonomic Design –

Three-Dimensional Custom-Design

Comfortable to wear without tight feeling

Three-Dimensional Custom-Design

Seamless fitting for all-round protection

Wider Elastic Nose Strip

Reduce the tenderness pain on the nose

Comfortable Ultra-Sonic Welding around The Chin

No irritation to the skin and comfortable to wear

Wider Knitted Ear Lanyard

Comfortable and breathable to extended wear



Innovative Material and Devices, Inc.

Address: 3F, Building #5, No. 615 Fengdeng Rd, Jiading District, Shanghai 201801, China

Toll free: +86 400 8202 169 Web: www.imdmedical.com